



ALUMINUM KIT

OCTA-SLIDE 2, 0-0.5 ppm

CODE 3569-01


QUANTITY	CONTENTS	CODE
50	*Aluminum #1 Tablets	*3943-H
50	*Aluminum #2 Tablets	*3944-H
2	Test Tubes, 5 & 10 mL, glass, w/cap	0898
1	Tablet Crusher	0175
1	Aluminum Octa-Slide 2 Bar, 0-0.5 ppm	7400-01
1	Octa-Slide 2 Viewer	1101

***WARNING:** Reagents marked with an * are considered to be potential health hazards. To view or print a Safety Data Sheet (SDS) for these reagents go to www.lamotte.com. Search for the four digit reagent code number listed on the reagent label, in the contents list or in the test procedures. Omit any letter that follows or precedes the four digit code number. For example, if the code is 4450WT-H, search 4450. To obtain a printed copy, contact LaMotte by email, phone or fax.

Emergency information for all LaMotte reagents is available from Chem-Tel: (US, 1-800-255-3924) (International, call collect, 813-248-0585).

To order individual reagents or test kit components, use the specified code number.

USE OF THE OCTA-SLIDE 2 VIEWER



The Octa-Slide 2 Viewer should be held so non-direct light enters through the back of the Viewer. Slide the Octa-Slide 2 Bar into the Viewer. Insert the reacted sample into the top of the Viewer. Match the color of the reaction to the color standards.

Warning! This set contains chemicals that may be harmful if misused. Read cautions on individual containers carefully. Not to be used by children except under adult supervision.










LaMOTTE COMPANY

Helping People Solve Analytical Challenges

PO Box 329 · Chestertown · Maryland · 21620 · USA
800-344-3100 · 410-778-3100 [Outside U.S.A.] · Fax 410-778-6394

Visit us on the web at www.lamotte.com

PROCEDURE

<p>1</p>  <p>Insert Aluminum Octa-Slide 2 Bar (7400-01) into the Octa-Slide 2 Viewer (1101)</p>	<p>2</p>  <p>Fill test tube (0898) to the 10mL line with sample water.</p>
<p>3</p>  <p>Add one *Aluminum #1 Tablet (3943). Crush with tablet crusher (0175) to disintegrate.</p>	<p>4</p>  <p>Add one *Aluminum #2 Tablet (3944). Crush with tablet crusher (0175) to disintegrate.</p>
<p>5</p> <p>Cap and mix.</p> 	<p>6</p> <p>Wait 10 minutes.</p> 
<p>7</p> <p>Mix the solution.</p> 	<p>8</p>  <p>Insert test tube into Octa-Slide 2 Viewer.</p>
<p>9</p>  <p>Match sample color to a color standard. Record as ppm Aluminum.</p>	