## Sulfid test paper



Sulfide test paper is an equal alternative to lead acetate paper, but with an important advantage, that it does not contain any hazardous materials, subject to declaration.

## Color reaction:

This white test paper changes, in the presence of sulfide ions and hydrogen sulfide, to brown-black, depending on the concentrations.

## Instructions for use:

Place a drop of the sample on the test paper, which will then change color to brown-black. Smaller quantities of sulfide now appear as a brown ring. For the determination of hydrogen sulfide, first wet the test paper with distilled water, and hold it up in the air to be tested. A brown coloration, on that part of the test

paper where water was applied, shows the presence of hydrogen sulfide.

## Limit of sensitivity:

The presence of 5 mg/L sulfide (S2-) produces an easily recognizable brown ring.

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany Tel.: +49 24 21 969-0 · Fax: +49 24 21 969-199 · info@mn-net.com · www.mn-net.com