

pH-Fix 2.0–9.0

en

Test strips for measuring pH values in the range of 2.0–9.0



IVD IVD for use by
qualified personnel only



Please refer to
instructions for use

REF 921 18

Pack content:

100 test strips
1 reference color scale pH 2.0–9.0

Application:

The pH-Fix test strips pH 2.0–9.0 are indicator strips for a variety of applications. The test strips are also suitable for measuring the pH value of gastric juice. Caution: The test strips pH-Fix 2.0–9.0 for measuring gastric juice should only be used by a qualified professional!

Please refer to a qualified professional for an interpretation of the test results.

General instructions for use:

Dip the test strip with all its test fields into the sample solution and read the pH value off the reference color scale. In weakly buffered solutions, keep the test strip immersed for a longer period of time (10–15 min).

Instructions for measuring the pH of gastric juice:

1. Collect sample in a clean vessel. Remove test strip and dip it with all its test fields into the sample until the color stops changing.
2. Remove the test strip and then remove any excess liquid by dabbing the strip along its edge on an absorbent cloth.
3. Compare the test fields with the reference color scale and read off the pH value (accuracy of reading: ± 0.2 pH).

Storage:

Avoid exposing the strips to sunlight and/or moisture. Keep container cool and dry (store below +30 °C).

Disposal:

Used test strips may be disposed of as domestic waste.



	Declaration of conformity (Product complies with the In Vitro Diagnostics Directive 98/79/EC of the European Union)	IVD	In Vitro Diagnostic Product
	Use by (not to be used after the use-by-date displayed on the packaging)	LOT	Batch identification
REF	Item number (for reordering purposes)		Not suitable for reuse. Dispose of after use.
	Storage temperature (4–30 °C) – Store at room temperature, and do not expose to direct sunlight.		Please refer to instructions for use! Prior to use, please read these instructions carefully.