# **REF 985091** 01.19 **Test 0-91** NANOCOLOR® Thiocyanate 50

en

Open test tube, add

Requisite accessories: piston pipette with tips

4.0 mL test sample (the pH value of the sample must be between pH 7 and 13), close and mix. Clean outside of test tube and measure immediately.

# Note:

Procedure:

This method can also be used for detecting a thiocvanate interference in the cvanide test. mg/L SCN<sup>-</sup> roughly corresponds to  $x \text{ ma/L } CN^-$ :

mg/L SCN⁻	mg/L CN⁻ (Test 0-31)
0.1	0.05
0.5	0.13
1.0	0.23
1.5	0.34
2.0	0.43

### Measurement:

For MACHEREY-NAGEL photometers see manual, test 0-91.

## Measurement when samples are colored or turbid:

For all NANOCOLOR® photometers see manual, use key for correction value.

#### Photometers of other manufacturers:

For other photometers check whether measurement of round glass tubes is possible. Verify factor for each type of instrument by measuring standard solutions.

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germanv

Method: Photometric determination as iron(III) thiocyanate

Range : Wavelength (HW = 5-12 nm): **Reaction time:** Reaction temperature:

#### Contents of reagent set:

20 test tubes Thiocyanate 50 1 test tube with blank value "NULL"

#### Hazard warning:

Test tubes contain Hydrochloric acid 10-25 %. For further information ask for a safety data sheet.

#### Interferences:

Nitrite, fluoride, anions of organic acids, phosphate, arsenate and borate interfere due to formation of complexes.

0.5-50.0 mg/L SCN<sup>-</sup>

470 nm

0 min 20-25 °C

The method can also be applied for the analysis of sea water after dilution (1+1).