

DIPYRIDYL Paper

en

for the rapid determination of ferrous ions

Colour reaction:

In the presence of Fe(II)-ions a red spot appears on white background.

Method of application:

Apply a drop of the weak mineral acid test solution (> pH 2) to the test paper. In the presence of Fe(II)-ions a red spot appears, or in the presence of small quantities of Fe(II)-ions a red ring.

Limit of sensitivity: 2 mg/L Fe(II)

Interferences:

The paper is specific for the determination of Fe(II)-ions. A few other metals result in pale discolourations of the test paper, but these are covered by the colour reaction of the ferrous ions. Fe(III)-ions result in a pale yellow reaction. In the presence of substantial quantities of Fe(III)-ions and small quantities of Fe(II)-ions, the Fe(III)-ion can be masked by the addition of small quantities of NaF.

MACHEREY-NAGEL GmbH & Co. KG

Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

Tel.: +49 24 21 969-0 · Fax: +49 24 21 969-199 · info@mn-net.com · www.mn-net.com