Test 1-42 01.19

NANOCOLOR® Fluoride



Method:

Photometric determination of fluoride with 1,8-dihydroxy-2-(4-sulfophenylazo)naphthalene-3,6-disulfonic acid (SPADNS)

 Cuvette:
 10 mm

 Range (mg/L F⁻):
 0.05–2.00

 Wavelength (HW = 5–12 nm):
 585 nm

 Reaction time:
 1 min (60 s)

 Reaction temperature:
 20–25 °C

Contents of reagent set:

8 x 75 mL Fluoride R1

Hazard warning:

Reagent R1 contains hydrochloric acid 10-25%.

For further information ask for a safety data sheet.

Preliminary tests:

If the order of magnitude of the concentration in a sample is not known, a preliminary test with Fluoride test paper 2–100 mg/L F⁻ (REF 90734) rapidly gives this information. From the order of magnitude the required dilution can be calculated and prepared directly.

Interferences:

Only repeatedly and thoroughly rinsed glassware should be used. Sea water and waste water require a distillation.

The following ions will not interfere: < 1000 mg/L Cu²+; < 500 mg/L Ca²+, Ni²+, Zn²+; < 200 mg/L Fe³+; < 100 mg/L SO₄²-; < 50 mg/L Cr(III); < 20 mg/L Si(IV); < 10 mg/L Cr(VI); < 5 mg/L PO₄³-, Cl₂; < 0,1 mg/L Al³+.

Procedure:

Requisite accessoires: volumetric flasks 25 mL, piston pipette with tips

Pour into two separate volumetric flasks:

Test sample	Blank value
20 mL test sample (the pH value of the sample	20 mL dist. water
must be between pH 7 and 13)	
3.0 mL R1, mix	3.0 mL R1, mix

Fill up sample and blank value to 25 mL mark with distilled water and mix again. After 1 min pour into cuvettes and measure.

Measurement:

For NANOCOLOR® photometers see manual, test 1-42.

Photometers of other manufacturers:

Verify calibration curve for each type of instrument by measuring standard solutions.

Analytical quality control:

NANOCONTROL Multistandard Metals 1 (REF 925015)

Analysis with reduced sample volume:

In order to increase the number of determinations, you can work with volumetric flasks of 10 mL: 8 mL test sample + 1.2 mL R1.

Disposal:

The contents of tubes and flasks can be washed into the drain with plenty of water.

Reference:

Standard methods for the examination of water and wastewater (4500-F⁻D)

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany
Tel.: +49 24 21 969-0 · Fax: +49 24 21 969-199 · info@mn-net.com · www.mn-net.com